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### 1. General description

Cables are intended to transmit electrical signals received from human organism to devices processing these signals - ECG apparatus.

Cables are non-sterile class I medical products.

In order to ensure proper cooperation of ECG cable with a monitor, choose an appropriate cable compatible with the device.

## **2. ECG cable characteristics**

The amount of conductors in the cable, cable length, method of connection of individual conductors to connectors, type of connectors, placement of electronic elements such as resistors and anti-surge elements in the cable, are adjusted to a specific type/model of apparatus, with which the cable is to cooperate.

Marking of the cable, method of its connection, purpose of individual tips and outlets, and other data, essential for correct conduction of examinations, is specified in the instruction of use of the device, with which the cable is intended to cooperate.

Cables are intended for use with undamaged patient's skin.

Cables should be used exclusively by qualified personnel.

## **3. Principles of use**

Before using the ECG cable, check the conformity of the product marking on the label and the product itself.

When using ECG cables, the following principles should be observed:

- cables should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the cable,
- place the cable so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin, clasp, latch ) when connecting and disconnecting the cable,
- protect the cable connectors against suffusion.

## **4. Installation**

Connect the cable pin to the appropriate ECG apparatus, turn on the device and verify the correctness of functioning according to monitor's instruction of use.

Pay attention to proper choice of the cable for an appropriate ECG apparatus.

## **5. Recommended cleaning, disinfection, and sterilisation**

These activities should be realized every time before use of the device and cable before the examination of the next patient. Before starting the cleaning, disinfection, or sterilisation, the cable should be disconnected from the device with which it cooperates.

### Cleaning

Wash the cable with warm soap water, delicately clean using soft cloth wetted with soap water (mild detergent). The cable should dry before use.

Do not use corrosive substances or other chemical preparations.

When cleaning, do not soak or immerse the cables, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to cable damage. Continuous bending of the conductors during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.

Avoid contact with strong solvents: aromatic, chloric, ketone, etheral, and ester. Do not allow for contact of cleaners with tips from the electrode side, and pins from ECG apparatus side.

### Disinfection

Disinfect cables by means of preparations dedicated to disinfect medical devices made of plastics based on active ingredients: isopropyl alcohol (70% concentration), glutaraldehyde, amines and quaternary ammonium compounds. The disinfection should be carried out in accordance with the instructions of the manufacturer of the preparation, and then leave the cable to dry. Do not allow for contact of preparat with tips from the electrode side, and pins from ECG apparatus side.

### Sterilization

The recommended method consists of sterilization using ethylene oxide. After sterilization and before use, cables should be accurately degassed in device intended for that purpose.

Sterilization with other methods is not recommended for it may permanently damage the cable.

## 6. Packaging and storing

ECG cables are packed individually. In order to maximize the product's working life, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:

- temperature: -25°C to + 55°C,
- relative humidity: 10% to 85%, do not allow for steam condensation
- protect against sunlight, moisture and dirt.

Cables in use should be stored in a dry place.

When not used, ECG cables should be loosely coiled, without sharp bends.

## 7. Environment conditions during work

- temperature +5°C to +45°C
- relative humidity 30% to 60%
- atmospheric pressure 700 + 1060 hPa

## 8. Warnings ⚠

Not observing the principles provided above may break electric circuit or lead to short-circuit. These signs are easy to identify by the ECG apparatus service staff. Damaged cable should be put away and replaced with an operational one. ECG cables should be used with fully functional devices, in the event of damage to the cable or device, stop the test immediately. Use only in rooms equipped with an efficient electrical and grounding installation.

ECG cables are intended for use with specified types of ECG apparatuses. The user is responsible for verification, before its use, of cable compatibility with the device. Incompatibility of these elements may lead to reduced accuracy and efficiency of measurements (or lack of them), and in extreme case - to damaging of the product. In order to ensure their compatibility, check the equipment cooperating with a given device in its instruction of use.

It is not recommended to use ECG cable in aroused patients, for movement negatively influences realized measurements. Relocation and positioning of the patient may cause disconnection of the cable from ECG electrodes and lack of readings.

Strong sources of electromagnetic field, such as those emitted by electrosurgical equipment, may negatively influence the correct functioning of the ECG cable. Do not use the ECG cable during magnetic resonance imaging (MRI).

Allergic reactions at contact of ECG cable with patient's undamaged skin are not observed. However, such reactions may not be excluded for persons particularly prone to allergies, the manufacturers are not known for this type of event. Improper cleaning or disinfection (or lack of them), and disregard for hygienic procedures may result in transfer of micro-organisms and cause infections.

When used during defibrillation, ECG cables should be used in conjunction with ECG equipment with defibrillation protection. For cables with defibrillation protection during the technical inspection, the effectiveness of defibrillation protection should be checked.

Do not use an ECG cable suspected of having mechanical or electrical defects, or packaging of which has been damaged. Suffused ECG cable should not be used. Do not change or modify the cable in any way. Such actions may influence correctness of functioning and accuracy of received readings.

Utilisation of used cable should be done in accordance with binding provisions.

## 9. Guarantee

A life time of the product is 10 years from the date of manufacture.

The guarantee period covers 1 year from the date of sale.

The general terms and conditions of the guarantee are described in a separate document available at the customer's request and at [www.sorimex.eu](http://www.sorimex.eu)

The guarantee does not include defects resulting from improper use and storage.

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- die EKG – Kabel für Diagnostik Typ W.KAB.KA; AW.KAB.KA; W.KAB.AD; AW.KAB.AD; W.KAB.AE; AW.KAB.AE; W.KAB.AP; AW.KAB.AP; W.KAB.AX; AW.KAB.AX; W.KAB.KB; AW.KAB.KB; W.KAB.KG; AW.KAB.KG; W.KAB.KS; AW.KAB.KS; W.KAB.KT; AW.KAB.KT
- die EKG - kabel für Überwachung Typ W.KAB.KA; AW.KAB.KA; W.KAB.AD; AW.KAB.AD; W.KAB.AE; AW.KAB.AE; W.KAB.AP; AW.KAB.AP; W.KAB.AX; AW.KAB.AX; W.KAB.KB; AW.KAB.KB; W.KAB.KM; AW.KAB.KM; W.KAB.AK; AW.KAB.AK; W.KAB.AS; AW.KAB.AS; W.KAB.Mm; AW.KAB.Mm; W.KAB.KS; AW.KAB.KS
- die Langzeit - EKG Kabel Typ W.KAB.KM; AW.KAB.KM; W.KAB.KH; AW.KAB.KH; W.KAB.KS; AW.KAB.KS; W.KAB.KB; AW.KAB.KB

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