SORIMEX

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EN

Instruction of use and storage disposable and reusable temperature sensors

1. General characteristics

Temperature sensors are designed for continuous measurement of patient temperature. This manual includes reusable and disposable temperature sensors. In terms of disposable sensors, this instruction applies with the exception of information on activities related to the re-use.

Temperature sensors are class IIb medical devices and are delivered non-sterile.

2. Characteristics of the sensor

The accuracy of the temperature measurement of the sensors is ± 0.1 °C for reusable sensors and ± 0.2 °C for single use sensors, in the range 25 °C + 45 °C.

The accuracy of the patient's vital function monitor is the sum of the accuracy of the temperature sensor and the accuracy of the device.

Compatibility - to ensure compliance and declared accuracy of equipment, SORIMEX temperature sensors should be used only with the specific equipment for which they were designed and labeled. Temperature sensors are interchangeable with other companies sensors, compatible with the same connector and matched temperature-resistivity characteristics.

Some sensor models end up on one side with a stainless, pointed disk tip - these are skin sensors.

To measure the skin temperature, apply the metal side of the temperature sensor.

3. Principles of use

- the temperature sensor and the cable must be cleaned before each use (not for single use items).
- in the case of prolonged use of the sensor on the patient's skin, the measurement site should be changed approximately every 4 hours to ensure the integrity of the skin of the patient,
- for patients who are excessively mobile, attach the sensors with a sticker.

4. Application

After making sure that the sensor is compatible with the device to be connected and that the plug of the sensor cable matches the socket on the device, insert the sensor plug into the device. Turn on the unit and check the correct operation according to the operating instructions of the device. Place the sensor tip in the measuring position and secure it. After starting the measurement, the readout will be performed after the readings have stabilized. Correct and safe use of the temperature sensor and its connecting cable requires systematic checks to be carried out at least once or several times a month, depending on the frequency of use and disinfection of the cord (not applicable to disposable sensors).

Check visually the insulation appearance, contact pins etc. Check the mechanical integrity of the connector. Do not use a sensor that appears to have mechanical or electrical defects and other signs of damage.

5. Suggestions for cleaning, disinfection, sterilization and instructions for use and storage

Disposable temperature sensors are designed for one patient, can not be reused and are not subject to cleaning, disinfection and sterilization procedures.

The following recommendations apply to reusable temperature sensors.

Cleaning

Cleanse the sensor and the connection cable with soap and warm water. Gently cleanse the sensor and the cable with the use of soft cloth moistened with water with soap or 70% isopropyl alcohol. Be careful not to damage the surface emitting the red and infrared radiation. The sensor and the cable should dry off before use. Do not use corrosive substances or chemicals other than 70% isopropyl alcohol. Do not soak or submerge the sensor. Excessive strength exercised during cleaning the sensor might cause the internal wires to tear. Those types of damage are not covered by the guarantee. Avoid any contact with strong dissolvents: aromatic, chlorinates, ketone, ether or ester. Do not allow the cleaners to have contact with the sensor's plug.

Disinfection

The sensor should be disinfected with 70 % isopropyl alcohol. Do not allow the alcohol to have contact with the plug from the plug side.

Sterilization

The preferred method of sterilization is with ethylene monoxide (procedure according to gas sterilizer's manufacturers). After sterilization and before use the sensors must be precisely and safely aired. The other methods of sterilization might permanently damage the sensor.

6. Packaging and storage

The temperature sensors are packed individually. The sensors should be stored in the original pack according to the storage conditions so that its operating time could be as long as possible.

Storage conditions:

- temperature: from 0°C to +40°C
- relative humidity: from 15% to 85%,
- protect from sunlight, dampness and dirt.

Keep the temperature sensors in a dry place.

The sensors and their cables, when not in use, should be loosely rolled up, without strong bends.

7. The conditions of the surroundings during operating

- temperature from +5°C to +45°C
- relative humidity from 30% to 75%
- atmospheric pressure 700 + 1060 hPa

8. Warnings ⚠

The temperature sensors are designed to be used with specific types of devices. Before the application, the user is responsible for checking the compatibility of the sensor and the device. The incompatibility of those elements might reduce the accuracy and efficacy of the measurement and in an extreme case damage the product. Check the manual of the camera, the equipment associated with it, in order to ensure their compatibility. Before starting to monitor the temperature, check the operation of the device to

which the sensor will be connected. Only used in conjunction with compatible devices ensure their compatibility.

Portable and mobile radiocommunication devices, electrosurgical devices, medical devices emitting electromagnetic waves can affect the operation and temperature indications.

Do not use a temperature sensor during MRI or CT tests. The conducted electricity might cause burns, and MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

Strong electromagnetic field emitted by the electro-surgical equipment might affect the proper functioning of the temperature sensors.

The temperature sensor is a product that has direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. The allergic reactions have not been observed when placing the sensor on the human skin, however such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people. Occurrence of allergic reactions or other skin reactions should be immediately subject to the physician.

Do not use a flooded temperature sensor in a defective package or one whose expiry date has been exceeded. Do not use the sensor if the sensor cable is damaged.

Remember to keep hygienic procedures when using the temperature sensor. If the skin is damaged during the skin preparation procedure for the temperature sensor application, select another application site.

To minimize the risk of skin irritation, replace the sensor with a new one on the skin.

Immediately after patient positioning and/or movement has been completed, check the temperature sensor connections with the temperature sensor cables and correct readings on the monitoring device.

Excessive movement of the patient adversely affects the measurements made, so disturbances in the reception of the signals are possible in the case of excitable patients. Temperature sensor cables should be arranged so as to eliminate accidental entanglement in the patient.

Temperature sensors can only be installed and repaired by trained and qualified personnel.

Utilization of the used temperature sensor should be in accordance with the applicable regulations.

8. Guarantee

The warranty period covers 1 year from the date of sale.

The shelf life of the disposable temperature sensor is 2 years from the date of manufacture.

This warranty does not cover defects resulting from improper use or storage.

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