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EN

Instruction of use and storage disposable SpO₂ sensors

1. General characteristics

SpO₂ sensors are used for non-invasive, constant measurement of the oxygenation of the arterial blood and the pulse.

There are various types of disposable SpO₂ sensors that differ in the structure of the products:

- SJR - the sensor placed on a foam tape with Velcro
- SJK - the sensor placed on a textile tape with adhesive layer
- SJP - the sensor placed on a foam tape with adhesive layer

The SpO₂ sensors mentioned in this manual are divided into four weight categories:

- adults - weight over 40 kilograms
- children - weight 10-50 kilograms
- infants - weight 3-20 kg
- neonates - weight under 3 kg

The SpO₂ sensors are active medical devices class IIb and are delivered in a non-sterile condition. They are type BF applied part and have IPX2 water protection.

The correct operating of the sensor with the device monitoring vital signs requires selecting a sensor compatible with the specific device. Compatibility is indicated on the product label.

2. Characteristics of the sensor

The accuracy of the saturation measurement:

±2% within the range 90% + 100%,

±3% within the range 70% + 90%

within the range under 69% - undefined.

The accuracy of the pulse measurement:

±2 within the 20 + 250 bpm range.

Compatibility:

To ensure accordance and declared accuracy of the products, the SpO₂ sensors should be used only with specific devices, for which they were designed and marked.

3. Principles of use

The SpO₂ sensor cannot be placed on the same arm as the SpO₂ cuff, arterial catheter or venous line. Before placing the product you have to remove the nail polish, since it can affect the accuracy of the measurement.

Avoid exposing the sensor to the intensive light source.

In case of hyperactive and overly mobile patients, adhesive sensors should be used or applied with the use of adhesive tape.

4. Application

1. Connect the SpO₂ sensor to the specific pulse oximeter, turn the device on and check if it operates correctly, according to the manual.
2. Choose the place of application on the patient's body. Preferred application spots:
 - adults and children: index finger, alternatively thumb, middle finger, ring finger or a big toe
 - infants and neonates: big toe, alternatively palm, wrist, ankle or foot
3. Place the sensor on a patient's body and position infrared connection. Fastening should be lasting but should not cause pressure.
4. Regularly monitor the skin condition in the place of the SpO₂ sensor's application to ensure correct blood circulation and integrity of patient's skin.

After the SpO₂ sensor has been placed you have to ensure the cable is aligned in a way that prevents the patient from entangling or breaking the product. The SpO₂ adapters that simplify the application of the sensor and cable's attachments such as bands, handles or clips are available in manufacturer offer.

The disposable SpO₂ sensor might be used by the same patient as long as the adhesive tape allows for the firm application of the product.

5. Transport and storage

The SpO₂ sensors are packed individually. The sensors should be stored in the original packaging accordingly to the storage conditions provided so the product does not get damaged during storage.

Transport and storage conditions:

- temperature: from 0°C to +40°C,
- relative humidity: from 15% to 85%, do not allow the steam to condense,
- pressure: 860 - 1060 hPa,
- protect from direct sunlight, dampness and dirt.

6. The conditions of the surroundings during operating

- temperature: from 0°C to +40°C
- relative humidity: from 15% to 85%
- atmospheric pressure: from 700 to 1060 hPa

7. Warnings

The SpO₂ sensors are designed to be used with specific types of devices. Before the application, the user is responsible for checking the compatibility of the sensor with the device. The incompatibility of those elements might reduce the accuracy and efficacy of the measurement and in an extreme case damage the product.

The accuracy of the sensor should not be checked with functional tester or oximeter's simulator. The manual of a given device and its equipment should be checked to ensure their compatibility. Only used in conjunction with compatible devices.

The use of the intravascular dye may distort the accuracy of the measurement.

The incorrect application of the SpO₂ sensor may distort the accuracy of the measurement as well.

The nail polish and/or tips should be removed before the application of the sensor since they may lead to inaccurate measurement.

The SpO₂ sensors should not be applied onto energized patient since their motion negatively affects measurement.

Do not use the SpO₂ sensor during MRI or CT tests. The conducted electricity might cause burns, and

MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

Mobile radio-communication devices might influence the readings of the pulse oximeter.

Strong electromagnetic field emitted by the electro-surgical equipment might affect the proper functioning of the SpO₂ sensors.

Strong sources of visible and infrared light such as direct sunlight, light from the devices measuring the level of bilirubin or infrared heat radiation might affect the functioning of the sensor and as a result distort the readings. In case either of the situations occur, cover the spot of the sensor's application with opaque, non-transparent material.

Do not use SpO₂ bands or adapters when placing the sensor or its cable on patient's body, which do not fit to the patient's anatomy. The additional pressure on the limb may result in alteration of venous pressure that may distort the findings.

The allergic reactions have not been observed when placing the sensor on the human skin, however such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people. Occurrence of allergic reactions or other skin reactions should be immediately subject to the physician.

In some cases there may be a need to regularly change the place of monitoring area, because the people skin's tolerance on the sensor might differ greatly.

When the monitor does not display the reliable pulse measurement it may be caused by improper place of the sensor's application or with deep pigmentation, which disables the normal light transmission. Then the place of the sensor's application must be changed or device dedicated for patient's with deep pigmentation of skin must be used.

The sensor that appears to have a mechanical or electrical defect, or its package is damaged must not be used. Flooded sensors must not be used as well.

It is not allowed to jerk/pull/stretch the SpO₂ sensor cable.

It is not allowed to make any modifications to the device. SpO₂ sensors can only be installed and repaired by trained and qualified personnel.

Utilization of the used SpO₂ sensor should be in accordance with the applicable regulations of electrical medical devices.

8. Garantie

The sensors shelf-life is two years.

The guarantee lasts for one year, beginning at the day of the purchase.

This warranty does not cover defects resulting from improper use or storage.

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The SpO₂ sensor generates, uses and may emit electromagnetic waves of radio frequency (RF). If the SpO₂ sensor is not used according to the directions specified in this instruction manual, it may cause electromagnetic interference. The user who uses the SpO₂ sensor should ensure that the sensor is used in an appropriate environment.

The SpO₂ sensors have been tested and recognized as compliant with restrictions specified in the EN 60601-1-2 standard for medical equipment. These restrictions can ensure suitable protection against electromagnetic interference (EMC), if the device is used in an environment which is suitable for the device and described in this instruction manual.

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Note 1: In case of frequencies of 80 MHz and 800 MHz, higher frequency range should be used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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a. It is not possible to theoretically calculate precisely the electromagnetic field intensities generated by stationary transmitters, such as base stations for radio telephony (mobile/wireless), by land-based movable transmitters, amateur transmitters, FM and AM radio broadcasting stations and TV broadcasting stations. To define a state of an electromagnetic environment resulting from the operation of the stationary radio frequency transmitters, it is necessary to make measurements on site. If measured electromagnetic field intensity at the location where the SpO2 sensor is used exceeds the admissible level, you should monitor whether the SpO2 sensor is operating properly. If improper operation of the SpO2 sensor is observed, further actions may be necessary, such as turning or moving the device.

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The SpO2 sensor is intended for using in electromagnetic environments characterized by controlled interference resulting from radio frequency emission. The buyer or user of the SpO2 sensor may prevent electromagnetic interference by maintaining the minimum distance between portable and movable communication devices using radio frequency (transmitters), and the SpO2 sensor, according to the following recommendations, depending on the maximum output power of the communication equipment.

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Rated output power of the transmitter (W) Maximale

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For transmitters characterized by maximum rated output power which has not been mentioned above, a recommended separating distance in meters (m) can be determined based on the equation suitable for the transmitter frequency, where P is a maximum rated output power of the transmitter in watts (W), specified by the manufacturer of the transmitter.

Note 1: For frequencies of 80 MHz and 800 MHz, the separating distance for higher frequency range is used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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