SORIMEX

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Suction electrodes

Saugelektrode

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EN

Instructions for use and storage of suction electrodes

These instructions are for reusable EP type ECG electrodes.

1. General characteristics

EP suction electrodes are intended to receive bioelectrical signals from the human body, in order to transfer them to the ECG device. Suction electrodes are reusable products.

2. Electrode characteristics

A suction electrode consists of:

- an electroconductive part, in contact with the patient's skin,
- a rubberbulb used to generate a vacuum retaining the electrode on the skin.

The suction electrode is connected to the ECG equipment using an ECG cable fitted with a Ø 4mm diameter banana plug, type H latch, or type M clamp.

3. Guidelines of use

When using ECG suction electrodes, the following guidelines principle should be followed:

- the skin at the point of electrode placement should not be covered with an ointment, lotion or other substances,
- in places of electrode's contact with patient's skin, the skin surface should be cleaned with a swob moistened with 70% ethyl alcohol,
- the connection cable should be arranged so as to prevent accidental movement of the electrode from the point of placement,
- in order to connect the electrode to ECG equipment, use cables compatible with a given type of ECG equipment,
- wires connected to the electrodes should be protected against accidental contact with other conductive parts, as well as earthed elements.

4. Application

- a) Choose the point of application of the limb electrodes on the patient's body. Remove lotion, gel, ointments, dirt and other non-conducting layers from the chosen points of application. If necessary, degrease the patient's skin using water and soap or other appropriate washing/disinfecting products, and then carefully dry the skin before placement of the electrode. Do not apply electrodes on moist, lubricated or oily skin, for this would make it impossible to establish a proper connection and signal transmission.
- b) Apply the gel intended for ECG examinations on the patient's skin at the point of contact with the electroconductive part of the ECG electrode, using enough to ensure contact with the whole metal surface of the electrode.
- c) Connect the ECG cable connectors to the suction electrodes.
- d) Place the electrodes on the patient's prepared skin.

5. Recommended cleaning, disinfection, and sterilization

The product is not intended to be sterilized by any available methods.

Cleaning and disinfection should be taken every time before using the electrodes for the examination of the next patient. Before starting the cleaning or disinfection, the suction electrodes should be disconnected from the ECG cable which they we connected to.

Cleaning

Wash the electrode with warm soapy water, delicately clean using a soft cloth. If necessary, use mild washing agents which will not react with the silver-covered electrode layer. The electrode should be dry before use. For washing, do not use agents containing abrasive materials, for example those intended to clean metals etc. When cleaning, do not soak or immerse the electrodes and do not wash under running water, do not use excessive strength or tools which may damage the electrode.

Disinfection

Suction electrodes should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry.

Warning: do not use ultrasonic washers.

6. Packaging and storing

Electrodes should be stored in their original packaging, in rooms at a temperature from 5°C to 30°C and air with relative humidity from 20% to 80%. Protect against sunlight, moisture and dirt.

Electrodes in use should be stored in conditions that have no negative effect on silver (in particular, in conditions that do not contain sulfur compounds).

7. Warnings

Electrodes are delivered to the client in a non-sterile state.

Before initiation of the ECG examination, it should be checked that the ECG device is functioning correctly. Use only together with compatible patient cables and ECG adapters.

Wires connected to the electrodes should be protected against accidental contact with other conductive parts, as well as earthed elements.

Electrodes should be connected to the ECG apparatus only by qualified medical personnel.

An ECG electrode is a product which has direct contact with the patient's skin; therefore allergic reactions may occur in patients prone to allergies. Occurrence of allergic reactions, or other skin-related reactions, requires immediate doctor's advice.

During the application of ECG electrodes, remember to apply sanitary procedures.

Do not apply electroconductive gel and electrodes on injured or irritated skin.

Immediately after positioning and/or moving the patient, check the connection of electrodes with the ECG cables/wires and the ECG device is showing the appropriate reading.

Do not use the electrode if it is wet or found in damaged packaging.

Excessive patient movement negatively influences the readings. Therefore, in the case of excitable patients, disruptions in ECG signals reception are possible.

Do not use ECG measurement electrodes immediately after defibrillation, for the time of signal recovery may be longer than 10 seconds.

Do not use ECG electrodes during Magnetic Resonance (MRI) examination or Computed Tomography (CT). Conducted electric current may inflict burns, and MRI and CT devices may cause inaccurate readings. An ECG electrode may influence MRI and CT images as well.

Disposal of used electrode should be done in accordance with legal requirements.

8. Shelf life

Shelf life depends on a proper continuity and quality of the layer of Ag / AgCl, which is covered with a sensor for contacting the skin of the patient. Improper cleaning can cause damage to the layer, which could result in signal

interference and inaccurate records ECG.

8. Garantie

The guarantee period covers 1 year from the date of sale.

The guarantee does not include defects resulting from improper use and storage and inadequate cleaning and disinfection.

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