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## Instruction of use and storage reusable SpO<sub>2</sub> sensors

### 1. General characteristics

SpO<sub>2</sub> sensors are used for non-invasive, constant measurement of the oxygenation of the arterial blood and the pulse.

There are various types of reusable SpO<sub>2</sub> sensors that differ in the structure of the products:

- SWA - adult/neonate universal wrap
- SWB - finger clip
- SWS - 'soft' sensor
- SWO - neonates band sensor
- SWY - multi-spot sensor
- SWU - ear clip sensor.

The SpO<sub>2</sub> sensors mentioned in this manual are divided into four weight categories:

- adults - weight over 40 kilograms
- children - weight 10-50 kilograms
- infants - weight 3-20 kg
- neonates - weight under 3 kg

The SpO<sub>2</sub> sensors are active medical devices class IIb and are delivered in a non-sterile condition. They are type BF applied part and have IPX2 water protection.

The correct operating of the sensor with the device monitoring vital signs requires selecting a sensor compatible with the specific device. Compatibility is indicated on the product label.

### 2. Characteristics of the sensor

The accuracy of the saturation measurement:

± 2% within the range 90% ÷ 100%,

± 3% within the range 70% ÷ 90%

within the range under 69% - undefined.

The accuracy of the pulse measurement:

±2 bpm within the 20 ÷ 250 bpm range.

#### Compatibility:

To ensure accordance and declared accuracy of the products, the SpO<sub>2</sub> sensors should be used only with specific devices, for which they were designed and marked.

### 3. Principles of use

The SpO<sub>2</sub> sensors and cables must be cleansed and disinfected before each and every use.

The SpO<sub>2</sub> sensor cannot be placed on the same arm as the SpO<sub>2</sub> cuff, arterial catheter or venous line. Before placing the product you have to remove the nail polish, since it can affect the accuracy of the measurement.

Avoid exposing the sensor to the intensive light source.

In case of hyperactive and overly mobile patients, adhesive sensors should be used or applied with the use of adhesive tape.

### 4. Application

1. Connect the SpO<sub>2</sub> sensor to the specific, compatible apparatus. Turn the device on and check if it operates correctly, according to the manual.
2. Choose the place of application on the patient's body. Preferred application spots:
  - adults and children: index finger, alternatively thumb, middle finger, ring finger, big toe or earlobe
  - infants and neonates: big toe, alternatively palm, wrist, ankle or foot

- Place the sensor on a patient's body and position infrared connection. Fastening should be lasting but should not cause pressure. For finger clip sensors, place the sensor on the finger, checking that the tip of the finger reaches the abutment part inside the sensor. Make sure the fingernail is under the abutment of the sensor clamp.
- Regularly monitor the skin condition in the place of the SpO<sub>2</sub> sensor's application.

After the SpO<sub>2</sub> sensor has been placed you have to ensure the cable is aligned in a way that prevents the patient from entangling or breaking the product. The SpO<sub>2</sub> adapters that simplify the application of the sensor and cable's attachments such as bands, handles or clips are available in manufacturer offer.

The safe and proper use of the SpO<sub>2</sub> sensor and its connection cable requires systematic checkups. The checkups must be conducted at least couple of times a month depending on the frequency of use and cable's disinfection.

Visually examine the front of the isolation, connection joints and pins etc. Examine the mechanical integrity of the joint. The sensor should not be used if there seems to be any mechanical or electrical fault

## 5. Suggestions for cleaning and disinfection

### Cleaning

Cleanse the sensor and the connection cable with soft cloth moistened with water with soap or mild detergent.

Be careful not to damage the surface emitting the red and infrared radiation. The sensor and the cable should dry of before use. Do not use corrosive substances or chemicals. Do not soak or immerse the sensor.

Excessive strength exercised during cleaning the sensor might cause sensor damage. Continuous bending of cables during use and cleaning might also lead the internal wires to tear. Those types of damage are not covered by the guarantee.

### Disinfection

The sensor should be disinfected with 70 % isopropyl alcohol. Disinfection should be done in accordance with the product manufacturer's instruction. Leave sensor to dry before use. Do not allow the alcohol to have contact with the plug from the joint side.

Avoid any contact with strong dissolvents: aromatic, chlorinates, ketone, ether or ester. Do not allow the cleaners to have contact with the sensor's plug.

## 6. Transport and packaging

The SpO<sub>2</sub> sensors are packed individually. The sensors should be stored in the original packaging accordingly to the storage conditions provided so the product does not get damaged during storage.

Transport and storage conditions:

- temperature: from -25°C to +55°C
- relative humidity: from 15% to 80%, do not allow the steam to condense
- pressure: 700 - 1060 hPa
- protect from direct sunlight, dampness and dirt

Keep the SpO<sub>2</sub> sensors in a dry place.

The sensors and their cables, when not in use, should be loosely rolled up, without strong bends.

## 7. The conditions of the surroundings during operating

- temperature from +0°C to +40°C
- relative humidity from 15% to 85%
- atmospheric pressure from 860 to 1060 hPa

## 8. Warnings

The SpO<sub>2</sub> sensors are designed to be used with specific types of devices. Before the application, the user is responsible for checking the compatibility of the sensor with the device. The incompatibility of those elements might reduce the accuracy and efficacy of the measurement and in an extreme case damage the product.

The accuracy of the sensor should not be checked with functional tester or oximeter's simulator. The manual of a given device and its equipment should be checked to ensure their compatibility. Only used in conjunction with compatible devices.

The use of the intravascular dye may distort the accuracy of the measurement.

The incorrect application of the SpO<sub>2</sub> sensor may distort the accuracy of the measurement as well.

The nail polish and/or tips should be removed before the application of the sensor since they may lead to inaccurate measurement.

The SpO<sub>2</sub> sensors should not be applied onto energized patient since their motion negatively affects measurement.

Do not use the SpO<sub>2</sub> sensor during MRI or CT tests. The conducted electricity might cause burns, and MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

Mobile radio-communication devices might influence the readings of the pulse oximeter.

Strong electromagnetic field emitted by the electro-surgical equipment might affect the proper functioning of the SpO<sub>2</sub> sensors.

Strong sources of visible and infrared light such as direct sunlight, light from the devices measuring the level of bilirubin or infrared heat radiation might affect the functioning of the sensor and as a result distort the readings. In case either of the situations occur, cover the spot of the sensor's application with opaque, non-transparent material.

Do not use SpO<sub>2</sub> bands or adapters when placing the sensor or its cable on patient's body, which do not fit to the patient's anatomy. The additional pressure on the limb may result in alteration of venous pressure that may distort the findings.

The allergic reactions have not been observed when placing the sensor on the human skin, however such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people. Occurrence of allergic reactions or other skin reactions should be immediately subject to the physician. In some cases there may be a need to regularly change the place of monitoring area, because the people skin's tolerance on the sensor might differ greatly.

When the monitor does not display the reliable pulse measurement it may be caused by improper place of the sensor's application or with deep pigmentation, which disables the normal light transmission. Then the place of the sensor's application must be changed or device dedicated for patient's with deep pigmentation of skin must be used.

The sensor that appears to have a mechanical or electrical defect, or its package is damaged must not be used. Flooded sensors must not be used as well.

For SWS type „soft“ finger sensors (soft silicone finger pad) it is not recommended to use excessive strength during application and when removing the sensor from the patient.

It is not allowed to jerk/pull/stretch the SpO<sub>2</sub> sensor cable.

It is not allowed to make any modifications to the device. SpO<sub>2</sub> sensors can only be installed and repaired by trained and qualified personnel.

Utilization of the used SpO<sub>2</sub> sensor should be in accordance with the applicable regulations of electrical medical devices.

## **9. Guarantee**

The warranty period covers 1 year from the date of sale.

Lifetime of the reusable sensor is indicated for 3 years.

This warranty does not cover defects resulting from improper use or storage.

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The SpO<sub>2</sub> sensor generates, uses and may emit electromagnetic waves of radio frequency (RF). If the SpO<sub>2</sub> sensor is not used according to the directions specified in this instruction manual, it may cause electromagnetic interference. The user who uses the SpO<sub>2</sub> sensor should ensure that the sensor is used in an appropriate environment.

The SpO<sub>2</sub> sensors have been tested and recognized as compliant with restrictions specified in the EN 60601-1-2 standard for medical equipment. These restrictions can ensure suitable protection against electromagnetic interference (EMC), if the device is used in an environment which is suitable for the device and described in this instruction manual.

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NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic waves are affected by absorption and reflection from structures, objects and people.

a. The ISM (industrial, scientific, and medical) bands from 150 kHz to 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands from 0.15 MHz to 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28, 0 MHz to 29.7 MHz, and 50.0 MHz to 54.0 MHz.

b. The field strengths of fixed transmitters, such as base stations of radio telephones (mobile/wireless) and terrestrial mobile radios, amateur radio, AM and FM radio stations and TV transmissions, cannot theoretically be predicted accurately. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the reusable SpO2 sensor is used exceeds the appropriate RF compliance level above, observe the reusable SpO2 sensor to verify normal operation. If abnormal operation is observed, additional measures such as changing the orientation or location of the reusable SpO2 sensor may be necessary.

c. In the frequency range of 150 kHz to 80 MHz, the field strength should be less than 3 V / m.

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For transmitters characterized by maximum rated output power which has not been mentioned above, a recommended separating distance in meters (m) can be determined based on the equation suitable for the transmitter frequency, where P is a maximum rated output power of the transmitter in watts (W), specified by the manufacturer of the transmitter.

Note 1: For frequencies of 80 MHz and 800 MHz, the separating distance for higher frequency range is used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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