SORIMEX

**Elektrody EKG**

**ECG Electrodes**

**EKG Elektroden**

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EN

## Instruction of use and storage of disposable ECG electrodes

Open the packaging immediately before usage.

The instruction is concerned with disposable, EK-S type ECG electrodes.

### 1. General characteristics

ECG electrodes are intended to receive bioelectrical signals from human organism, in order to transfer them to the ECG device.

### 2. Electrode characteristics

After removal of the foil cover from, the electrode is placed on the undamaged skin at the spot from which the electrical signal is to be received.

A terminal used to connect the electrode to the ECG device is located on the external side of the electrode.

Parts of electrodes connecting to the patient's skin (adhesive layer and gel) are made of materials neutral for human organism. Electrodes marked as RTG are radiotranslucent in X-rays.

### 3. Principles of use

When using ECG electrodes, the following principles should be followed:

- The skin at the spot of electrode placement should not be covered with an ointment, lotion etc. It is recommended to remove hair from these spots.
- The connection cable should be arranged so as to prevent accidental detachment of the electrode.
- The electrode should be connected to ECG devices with cables designed to be compatible with given type of ECG device.
- In some types of electrodes, the connection cable is supposed to be inserted in an oblong opening or incision on the electrode. Details regarding the method of application of such electrodes are described in the point 4.1. of this instruction.
- Connection cables connected to the electrodes should be protected from accidental contact with other current-conducting parts, as well as with earthed elements.
- After the examination is finished, electrodes should be disposed of in a way appropriate for medical products.
- If electrodes are used over several days, firm placement of electrodes and absence of skin irritation should be confirmed on a daily basis. Once applied, the location of electrodes must not be changed.

### 4. Application

#### 4.1. Application of standard ECG electrodes

- a) Open (slit) the pouch along the line marked on the packaging.
- b) Choose the place of application on patient's body. Remove lotion, gel,


ointments, dirt and other non-conducting layers from chosen places of application. If needed, clean patient's skin using water and soap or other appropriate detergent or sanitizer, and then carefully dry up the skin before placement. Do not stick to moist, lubricated or oily skin, for it would make it impossible to establish a stable connection and signal transmission. Before application of the electrode, always confirm that the place of application is dry. Do not stick to overly hairy spots.

- c) Connect the drop cable with the electrode placed on the plastic cover.
- d) After removing the transparent, plastic cover, apply electrodes on previously prepared skin. From the sticky side of the electrode unglue the white paper which facilitates removing electrode from the transparent foil. If needed, smooth the surface of the electrode with finger, gently pressing the edges.
- e) After the examination is finished, the electrode should be levered in order to gently unglue it from the patient's skin.

#### 4.2. Application of Holter examination ECG electrodes

Points from a) to c) proceed as described above.

- d) When preparing electrodes for application, you should ensure that the in-cut for the cable is facing **DOWNWARDS**.
  - e) Without ungluing ECG electrodes from the plastic cover, pull the cable through the cut opening - according to the picture - from the glue side to the external side of the electrode.
- 
- f) Connect the drop cable placed in the slit with the electrode. Fastener/clasp lead out should be pointed upwards. Leave about 3-4 centimetres of loose cable (above the opening) in order to prevent its disconnection from the electrode.
  - g) After ungluing from the plastic cover, place the electrode on a dry and previously cleaned skin.
  - h) Unglue the safety belt located below the electrode opening and press the electrode as to immobilize the previously prepared connection cable.
  - i) If needed, smooth the surface of the electrode with finger, gently pressing the edges.
  - j) After the examination is finished, the electrode should be levered in order to gently unglue it from the patient's skin.

#### 5. Packaging and storing

Electrodes should be stored in the original packaging, in temperature ranging from 5°C to 30°C, protected from the sunlight, dampening and dirt.

After the packaging is opened, unused electrodes should be left in the original packaging and secured by folding an edge of the packaging. Opened can be stored for max. 5 - 7 days.

#### 6. Warnings

Before initiation of the ECG examination, the correct functioning of ECG device should be ensured. Use only together with compatible patient's cables and ECG adapters. For disposable ECG electrodes, appropriate types of cable endings consist of fasteners and clasps. In the case of cables ending with a banana plug, an appropriate ECG adapter should be used.

Do not use ECG electrodes during Magnetic Resonance (MRI) examination or Computed Tomography (CT). Conducted electric current may inflict burns, and MRI and CT devices may cause inaccurate readings.<sup>6</sup> EKG electrode may influence MRI and CT images as well.

Computed Tomography (CT). Conducted electric current may inflict burns, and MRI and CT devices may cause inaccurate readings. EKG electrode may influence MRI and CT images as well.

Strong sources of electromagnetic field, such as those emitted by electrosurgical equipment, may negatively influence the correct functioning of the ECG electrode. In the case of disposable ECG electrode usage in operating theatres, cables and equipment with appropriate protection against electromagnetic disruptions should be used.

Do not use the electrode if it is suffused, found in damaged packaging, or which expiry date has passed.

Gel of electrodes stored in open packaging may dry up, so the pouch should not be opened if there is no necessity of conducting the examination. Electrodes from opened packaging should be used as soon as possible.

An ECG electrode is a product which has a direct contact with patient's skin; therefore allergic reactions may occur in patients prone to allergies. If such allergic reactions, or other skin-related reactions, occur, they should be immediately examined by a doctor.

In view of firmly adhesive glue, the product is not appropriate to be used in patients with sensitive skin. In order to minimize the risk of skin irritation, when changing an electrode for a new one, the product should be placed on different spot on the skin.

Caution is advised in the use of abrasive products (pastes, gels, creams) preparing the patient's skin to ECG examination, especially in patients with hypersensitive skin.

During the preparation of ECG electrodes, sanitary procedures should be followed. Recommended method of patient's skin preparation before application of an electrode: shave overly hairy skin, clear up with several gentle movements with gauze or swab impregnated with a detergent or sanitizer, do not injure the epidermis, dry up. In the case of damaging of the epidermis during the preparation for electrode application procedure, other place of application should be chosen. Do not use alcohol solution or petrol solution to prepare the skin before ECG electrodes application.

Immediately after positioning and/or relocating the patient, both connection of electrodes to EKG cables and accuracy of the ECG device reading should be checked.

Excessive patient movement negatively influences the readings. Therefore, in the case of excitable patients, disruptions in ECG signals reception are possible.

The product is not intended for internal use. In the case of accidental swallowing of the electrode, a doctor should be consulted immediately.

Do not allow for contact with external electric current!

Utilisation of used ECG electrode should be done in accordance with regulations concerning medical products.

## **7. Shelf life**

The shelf life is 2 years from manufacture date under recommended storage conditions (section 5) in closed original packaging.

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