SORIMEX

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Instruction of use and storage disposable and reusable temperature sensors

EN

1. General characteristics

Temperature sensors are designed for continuous measurement of patients body temperature. This manual includes reusable and disposable temperature sensors. In terms of **disposable sensors**, this instruction applies **with the exception of information on activities related to the re-use**.

There are various types of temperature sensors that differ in the structure of the products:

- skin sensor for adults,
- skin sensor for children,

- general sensor for adults,
- general sensors for children.

The temperature sensors mentioned in this manual are divided into following weight categories:

- adults - weight over 40 kilograms
- children - weight 10-50 kilograms

The temperature sensors are active medical devices class IIb and are delivered in a non-sterile condition. They are type BF applied part and have IPX2 water protection. Intended for short-term use with intact patient skin.

Reusable temperature sensor models with the letter "O" as the last character of their reference number are invasive medical devices for use in natural orifices and these with letter „S" - for measuring temperature from the surface of intact patient's skin.

The correct operating of the sensor with the device (monitor) requires selecting a sensor compatible with the specific device. Compatibility is indicated on the product label.

2. Characteristics of the sensor

The accuracy of the temperature measurement of the sensors is:

- ± 0.2 °C for single use sensors,
- ± 0.1 °C for reusable sensors,

in the range 25 °C + 45 °C.

The accuracy of the patient's vital function monitor is the sum of the accuracy of the temperature sensor and the accuracy of the device.

Compatibility

To ensure compliance and declared accuracy of equipment, SORIMEX temperature sensors should be used only with the specific equipment for which they were designed and marked on product label. Temperature sensors are interchangeable with other companies sensors, compatible with the same connector and matched temperature-resistivity characteristics.

3. Principles of use

- Some sensor models are ended up on one side with a stainless, pointed disk tip - these are skin sensors. To measure the temperature, apply the metal side of the temperature sensor to the patient's skin.
- The temperature sensor and the cable must be cleaned before each use (not for single use items).
- In the case of prolonged use of the sensor (disposable as well as reusable) on the patient's skin, the measurement site should be changed approximately every 4 hours to ensure the integrity of the patient's skin.
- Disposable temperature sensors should not be re-used for other patients.
- For patients who are excessively mobile, attach the sensors with a sticker.

4. Application

After ensuring that the sensor is compatible with the device to be connected and that the plug of the sensor cable matches the socket on the device, insert the sensor plug into the device. Turn on the unit and check the correct operation according to the operating instructions of the device. Place the sensor tip in the measuring position and secure it. After starting the measurement, the readout should be performed after the readings have stabilized. Correct and safe use of the temperature sensor and its connecting cable requires systematic checks to be carried out. They should be conducted at least once or several times a month, depending on the frequency of use and disinfection of the cord (not applicable to disposable sensors).

5. Suggestions for cleaning and disinfection

Disposable temperature sensors are designed for one patient. They can not be reused and are not subjected to cleaning and disinfection procedures.

The following recommendations apply to reusable temperature sensors.

Cleaning

Gently cleanse the sensor and the connection cable with the use of soft cloth moistened with water and soap or delicate detergent. Be careful not to damage the surface of the sensor. The sensor and the cable should dry of before use. Do not use corrosive or other aggressive chemicals substances. Do not soak or submerge the sensor.

Excessive strength exercised during cleaning the sensor might cause the internal wires to tear, leading to sensor damage. Constant bending of cables during cleaning and usage may also lead to internal wires tear. Those types of damages are not covered by the guarantee.

Disinfection

Sensors should be disinfected with 70% isopropyl alcohol and left to dry out.

Avoid any contact with strong dissolvents: aromatic, chlorinated, ketone, ether or ester. Do not allow the cleaners and disinfectants to have contact with the sensor's plug.

6. Transport and storage

The temperature sensors are packed individually. The sensors should be stored in the original packaging accordingly to the storage conditions provided so the product does not get damaged during storage.

Transport and storage conditions:

- temperature: from 0°C to +40°C
- relative humidity: from 15% to 85%,
- protect from direct sunlight, dampness and dirt.

Keep the temperature sensors in a dry place.

The sensors and their cables, when not in use, should be loosely rolled up, without strong bends.

7. The conditions of the surroundings during operating

- temperature: from +15°C to +35°C,
- relative humidity: from 15% to 85%.

8. Warnings

The temperature sensors are designed to be used with specific types of devices. Before the application, the user is responsible for checking the compatibility of the sensor with the device. The incompatibility of those elements might reduce the accuracy and efficacy of the measurements. Check the manual of the monitor, the equipment associated with it, in order to ensure their compatibility. Before starting to monitor the temperature, check the operation of the device to which the sensor will be connected. Only used in conjunction with compatible devices.

Portable and mobile radiocommunication devices, electrosurgical devices, medical devices emitting electromagnetic waves can affect the operation and temperature indications.

Do not use a temperature sensor during MRI or CT tests. The conducted electricity might cause burns, and MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

The temperature sensor is a product that has direct contact with the patient's skin, and therefore allergic reactions may occur in especially sensitive individuals. The allergic reactions have not been observed when placing the sensor on the human skin. However, such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people. Occurrence of allergic reactions or other skin reactions should be immediately subject to the physician.

Do not use a flooded or in a defective package temperature sensor. Disposable sensors should not be used after expiry date has been exceeded. Do not use the sensor if the sensor cable is damaged.

Check visually the insulation appearance, contact pins etc. Check the mechanical integrity of the connector. Do not use a sensor that appears to have mechanical or electrical defects and other signs of damage.

Remember to keep hygienic procedures when using the temperature sensor. If the skin is damaged during the skin preparation procedure for the temperature sensor application, select another application site.

To minimize the risk of skin irritation, when replacing the sensor with a new one, select another

application place on the patient's skin.

Immediately after patient positioning and/or movement has been completed, check the temperature sensor connections with the temperature sensor cables and correct readings on the monitoring device.

Excessive movement of the patient adversely affects the measurements made, so disturbances in the reception of the signals are possible in the case of excitable patients.

Temperature sensor cables should be arranged so that eliminates accidental entanglement of the patient.

It is not allowed to make any modifications to the device. Temperature sensors can only be installed and repaired by trained and qualified personnel.

Utilization of the used temperature sensor should be in accordance with the applicable regulations for electrical medical devices.

9. Guarantee

The warranty period covers 1 year from the date of sale.

Use disposable temperature sensor by 2 years from production date.

Lifetime of the reusable temperature sensor is 3 years.

This warranty does not cover defects resulting from improper use or storage.

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The temperature sensor is intended for use in an electromagnetic environment which is described below. The buyer or user of the temperature sensor should ensure using the device in such environment.

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| | | | |
|-------------------------------|--------------------------------|----------|--|
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2,7 GHz | Complies | Electromagnetic field intensity off the screen from stationary radio transmitters measured during electromagnetic field tests on site, should not exceed 3 V/m. Interference may be available near devices marked with:  |
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | Complies | |

The method of calculating the distance is given below. If a known transmitter is present, the specified distance can be calculated using equations.

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Note 1: In case of frequencies of 80 MHz and 800 MHz, higher frequency range should be used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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For transmitters characterized by maximum rated output power which has not been mentioned above, a recommended separating distance in meters (m) can be determined based on the equation suitable for the transmitter frequency, where P is a maximum rated output power of the transmitter in watts (W).

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