SORIMEX

PL Kable SpO2

EN SpO2 cables

DE SpO2-Kabel

SORIMEX sp. z o.o. (dawniej SORIMEX sp. z o.o. sp. k.)
Równinna 25, 87-100 Toruń, Poland
office@sorimex.eu ♦ www.sorimex.eu

Rev. 2025-01-07

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PL

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- **kabli SpO2**

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This instruction for use concerns:

- **SpO2 cables**

1. Intended use and operation of the product

An SpO2 cable is intended to be used for connecting a saturation sensor with a device used to measure the blood oxygen saturation.

The cable is used to transmit supply voltage from the device to the sensor and voltage from the sensor that is used by the device to determine the blood oxygen saturation and the heart rate. The cable consists of a multi-wire cable of an appropriate length ending with connectors suitable for the saturation sensor and the device with which it is to work.

2. Target patient groups

SpO2 cables are intended to be used by qualified medical personnel in neonates, infants, children and adults.

3. Contraindications

No known contraindications to the use of the product.

4. Possible adverse effects

SpO2 cables can contact the uninjured patient's skin. No allergic reactions are observed in contact with the patient's skin; however, such reactions cannot be completely excluded in patients particularly susceptible to allergies.

5. Warnings

- 5.1.** A damaged SpO2 cable should be disconnected and replaced with one in a good operating condition. SpO2 cables must be used with units that are in a perfect operating condition.
- 5.2.** SpO2 cables are intended to be used with specific types of saturation sensors and devices. Before each use of the cable, check whether a given cable is suitable for a given type of the device and the sensor, with which it is to operate.
- 5.3.** When installing the SpO2 cable, follow the operating manual of the device.

6. Safety precautions

- 6.1.** Do not use an SpO2 cable with any mechanical or electrical defects.
- 6.2.** Do not use SpO2 cables that were in contact with any liquid.
- 6.3.** SpO2 cables should not be bent excessively (the bend radius should be at least 60 mm).
- 6.4.** Do not pull or yank a cable.
- 6.5.** A cable should be routed in such way that there is no risk of its mechanical damage caused, e.g., by stepping on, driving over, placing a heavy/sharp object on, or cutting it.

- 6.6. When connecting or disconnecting a cable, hold it by a connection (socket, plug) housing.
- 6.7. Movement, transfer and positioning of the patient may interfere with reception of collected signals and affect measurement results.
- 6.8. Do not use the SpO₂ cable during Magnetic Resonance Imaging (MRI) or Computed Tomography (CT) scans.

7. How to use the cable

- 7.1. Each time before using the product, clean and disinfect it.
- 7.2. Select an SpO₂ cable appropriate for the sensor and the device with which it is to operate.
- 7.3. Connect the SpO₂ cable plug to the device socket, and on the other side, connect the SpO₂ cable socket to a plug of the saturation sensor.
- 7.4. Switch the device on and check it for correct operation in accordance with the monitor operating manual.

8. Cleaning and disinfection

These actions should be performed each time before the SpO₂ cable is used. Before cleaning or disinfection, a cable should be disconnected from a device and a sensor with which it operates.

8.1. Cleaning

A cable or a connecting wire should be cleaned by delicately wiping it with a soft cloth wetted with warm water and soap or mild detergent. Before the use, the cable/wire should be dry.

Do not use any corrosive agents or other aggressive chemical formulations. Avoid any contact with strong solvents: aromatic, chlorine, ketone, ethers or esters.

Do not allow any contact of cleaning agents with a terminal of a cable on the sensor side and with plugs on the device side.

8.2. Disinfection

Disinfect cables by wiping them delicately with formulations dedicated to disinfect plastic medical devices, based on active ingredients such as (70%) isopropyl alcohol, glutaraldehyde, amines, and quaternary ammonium compounds. Disinfect according to instructions of a formulation manufacturer, and leave the cable to dry.

9. Signs of product degradation

Before each use, perform visual inspection of the product. SpO₂ cables cannot be used any longer when any of the below is noticed:

- breaks in the cable jacket (external insulation) continuity, revealing conducting elements of the cable (wires, shield, connecting layer)
- deep scratches not affecting the jacket continuity, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection
- cracks in plastic material covering product components, making electrical connections visible

- deep scratches in plastic material covering product components that do not make electrical connections visible, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection of the product
- partial or complete tearing of cables from connections/plugs or breaking of the jacket in those areas

Any of the above defects disqualifies the product from further use. A product with such defects should either be repaired or disposed of.

10. Storage conditions

10.1. Storage conditions:

- temperature: -25°C to +55°C
- relative humidity: 15% to 85%, avoid condensation
- atmospheric pressure: 700 to 1060 hPa
- protect against sunlight, moisture and contamination

10.2. SpO₂ cables should be stored in the specified storage conditions, to prolong the product life as far as possible.

10.3. During storage, the SpO₂ cables should be loosely wound, without kinks.

11. Compatibility

SpO₂ cables are intended as intermediary elements for connections between SpO₂ sensors and devices.

The indication of the device operating in conjunction with the product is given on the label unit packaging of the product and on the product itself

12. Guarantee

SpO₂ cables are intended to be used as intermediate components connecting SpO₂ sensors with devices.

A device operating together with the product is specified in a label on the product primary packaging and on the product itself.

13. Description of symbols used

Product catalogue number

Pieces per packaging

Production batch number

Medical device

Production date

Manufacturer

Consult the instruction for use

Caution

Protect against sunlight

CE mark (conformance to the Regulation EU 2017/745)

Keep dry

Acceptable temperature limits

Non-sterile device

Atmospheric pressure limitation

Reduced moisture

14. Notification of incidents related to the product use

Each serious incident related to the product should be notified to the manufacturer and to an authority of the Member State competent for the user or patient's place of residence.

15. Product disposal

The used cable should be disposed of in accordance with internal hygiene procedures in force at a given health care facility and with local regulations.

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