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EN

Instruction of use and storage disposable SpO2 sensors

1. General characteristics

SpO2 Sensors are used for constant, approximate and non-invasive measurement of the oxygenation of the arterial blood and the pulse.

There are various types of disposable SpO2 sensors that differ in the structure of the products:

- SJR - the sensor placed on a foam tape with Velcro
- SJK - the sensor placed on a textile tape with adhesive layer
- SJP - the sensor placed on a foam tape with adhesive layer.

The SpO2 sensors mentioned in this manual are divided into four weight categories:

- adults - weight over 30 kilograms
- children - weight 15-40 kilograms
- infants - weight 1-15 kg
- neonates - weight under 3 kg

The correct operating of the sensor with the device (pulse oximeter) requires selecting a sensor compatible with the specific device.

The disposable SpO2 sensors belong to the IIb class and are delivered in a sterile condition.

2. Characteristics of the sensor

The accuracy of the saturation measurement:

- adults: $\pm 2\%$ within the range $70\% \div 100\%$, undefined within the range under 70%
- neonates: $\pm 3\%$ within the range $70\% \div 100\%$, undefined within the range under 70%

The accuracy of the pulse measurement: ± 3 within the $20 \div 250$ bpm range.

Compatibility: to ensure accordance and declared accuracy of the products, the SpO2 sensors should be used only with specific devices, for which they were designed and marked.

Mechanical integrity: the SpO2 sensors are made from the best quality materials so that they withstand even severe operating conditions. The SpO2 sensors do not have any uncovered metal parts.

3. Principles of use

The SpO2 sensor cannot be placed on the same arm as the SpO2 cuff, arterial catheter or venous line. Before placing the product you have to remove the nail polish, since it can affect the accuracy of the measurement. The sensor cannot be used inside or near the MRI. Avoid exposing the sensor to the intensive light source. In case of the long use, the place of application should be changed once every 4 hours to ensure the integrity of the patient's skin. In case of hyperactive and overly mobile patients, adhesive sensors should be used or applied with the use of plaster.

4. Application

1. Connect the SpO2 sensor to the specific pulse oximeter, turn the device on and check if it operates correctly, according to the manual.
2. Choose the place of application on the patient's body. Preferred application spots:
 - adults and children: index finger, alternatively thumb, middle finger, ring finger or a big toe.
 - infants and neonates: big toe, alternatively palm, wrist, ankle or foot.
3. Place the sensor on a patient's body and position infrared connection. Fastening should be lasting but should not cause pressure.
4. Regularly monitor the skin condition in the place of the SpO2 sensor's application to ensure correct blood circulation and integrity of patient's skin.

After the SpO2 sensor has been placed you have to ensure the cable is aligned in a way that prevents the patient from entangling or breaking the product. The SpO2 adapters that simplify the application of the sensor and cable's attachments such as bands, handles or clips are available on stock.

The disposable SpO2 sensor might be used by the same patient as long as the adhesive tape allows for the firm application of the product.

5. Packaging and storage

The SpO2 sensors are packed individually. The sensors should be stored in the original pack according to the storage conditions so that its operating time could be as long as possible.

Storage conditions:

- temperature: from 0°C to +50°C
- relative humidity: from 0% to 60%, do not allow the steam to condense

6. The conditions of the surroundings during operating

- temperature from +10°C to +40°C
- relative humidity from 0% to 60%
- atmospheric pressure 700 ÷ 1060 hPa

7. Warnings

The SpO2 sensors are designed to be used with specific types of pulse oximeters. Before the application, the user is responsible for checking the compatibility of the sensor and the pulse oximeter. The incompatibility of those elements might reduce the accuracy and efficacy of the measurement and in an extreme case damage the product. The accuracy of the sensor should not be checked with functional tester or oximeter's simulator. The manual of a given device and its equipment should be checked to ensure their compatibility.

The use of the intravascular dye may distort the accuracy of the measurement. The incorrect application of the SpO2 sensor may distort the accuracy of the measurement as well. The nail polish and/ or tips should be removed before the application of the sensor since they may lead to inaccurate measurement. The SpO2 sensors should not be applied onto energized patient since their motion negatively affects measurement.

Do not use the SpO2 sensor during MRI or CT tests. The conducted electricity might cause burns, and MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

Strong electromagnetic field emitted by the electro-surgical equipment might affect the proper functioning of the SpO2 sensors.

Strong sources of visible and infrared light such as direct sunlight, light from the devices measuring the level of bilirubin or infrared heat radiation might affect the functioning of the sensor and as a result distort the readings. In case either of the situations occur, cover the spot of the sensor's application with opaque, non-transparent material.

Do not use other SpO2 bands or adapters than indicated by the manufacturer when placing the sensor or its cable on patient's body.

The additional pressure on the limb may result in alteration of venous pressure that may distort the findings.

The allergic reactions have not been observed when placing the sensor on the human skin, however

such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people. In some cases there may be a need to regularly change the place of monitoring area, because the people skin's tolerance on the sensor might differ greatly.

Mobile radio-communication devices might influence the readings of the pulse oximeter. When the monitor does not display the reliable pulse measurement it may be caused by improper place of the sensor's application, or the monitored area is too thick, too thin or with deep pigmentation, which disables the normal light transmission. Than the place of the sensor's application must be changed.

The sensor that appears to have a mechanical or electrical defect, or its package is damaged must not be used.

Inundated sensors must not be used as well.

Do not alter or modify the sensor in any way. Such actions may influence the preceisness of readings.

The used sensor should be disposed of in accordance with local requirements.

8. Garantie

The sensors shelf-life is two years.

The guarantee lasts for one year, beginning at the day of the purchase.

Guarantee does not cover the improper use and modification of the sensors.

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The SpO₂ sensor generates, uses and may emit electromagnetic waves of radio frequency (RF). If the SpO₂ sensor is not used according to the directions specified in this instruction manual, it may cause electromagnetic interference. The user who uses the SpO₂ sensor should ensure that the sensor is used in an appropriate environment.

The SpO₂ sensors have been tested and recognized as compliant with restrictions specified in the EN 60601-1-2 standard for medical equipment. These restrictions can ensure suitable protection against electromagnetic interference (EMC), if the device is used in an environment which is suitable for the device and described in this instruction manual.

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The SpO2 sensor is intended for use in an electromagnetic environment which is described below. The buyer or user of the SpO2 sensor should ensure using the device in such environment.

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Manufacturer's recommendations and statement – electromagnetic resistance

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a. It is not possible to theoretically calculate precisely the electromagnetic field intensities generated by stationary transmitters, such as base stations for radio telephony (mobile/wireless), by land-based movable transmitters, amateur transmitters, FM and AM radio broadcasting stations and TV broadcasting stations. To define a state of an electromagnetic environment resulting from the operation of the stationary radio frequency transmitters, it is necessary to make measurements on site. If measured electromagnetic field intensity at the location where the SpO2 sensor is used exceeds the admissible level, you should monitor whether the SpO2 sensor is operating properly. If improper operation of the SpO2 sensor is observed, further actions may be necessary, such as turning or moving the device.

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Recommended separation distances between portable and cellular communication devices and the SpO2 sensor

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The SpO2 sensor is intended for using in electromagnetic environments characterized by controlled interference resulting from radio frequency emission. The buyer or user of the SpO2 sensor may prevent electromagnetic interference by maintaining the minimum distance between portable and movable communication devices using radio frequency (transmitters), and the SpO2 sensor, according to the following recommendations, depending on the maximum output power of the communication equipment.

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