## Czujniki SpO<sub>2</sub>

SpO<sub>2</sub> sensors

SpO<sub>2</sub> Sensoren

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**EN**

## **Instruction of use and storage reusable SpO2 sensors**

### **1. General characteristics**

SpO2 Sensors are used for constant, approximate and non-invasive measurement of the oxygenation of the arterial blood and the pulse.

There are various types of disposable SpO2 sensors that differ in the structure of the products:

- SWB - finger clip for adults
- SWS - 'soft' sensor
- SWO - band sensor
- SWY - multi-spot sensor
- SWA - alternative wrap
- SWU - ear clip.

The SpO<sub>2</sub> sensors mentioned in this manual are divided into four weight categories:

- adults - weight over 30 kilograms
- children - weight 15-40 kilograms
- infants - weight 1-15 kg
- neonates - weight under 3 kg

The correct operating of the sensor with the device (pulse oximeter) requires selecting a sensor compatible with the specific device.

The reusable SpO<sub>2</sub> sensors belong to the IIb class and are delivered in a sterile condition.

## 2. Characteristics of the sensor

The accuracy of the saturation measurement:  $\pm 3\%$  within the range  $70\% + 100\%$ , undefined within the range under  $70\%$ .

The accuracy of the pulse measurement:  $\pm 3$  within the  $20 + 250$  bpm range.

Compatibility: to ensure accordance and declared accuracy of the products, the SpO<sub>2</sub> sensors should be used only with specific devices, for which they were designed and marked.

Mechanical integrity: the SpO<sub>2</sub> sensors are made from the best quality materials so that they withstand even severe operating conditions. The SpO<sub>2</sub> sensors do not have any uncovered metal parts.

Safety – type of electrocution protection – BF

## 3. Principles of use

The SpO<sub>2</sub> sensors and cables must be cleansed and disinfected before each and every use.

The SpO<sub>2</sub> sensor cannot be placed on the same arm as the SpO<sub>2</sub> cuff, arterial catheter or venous line. Before placing the product you have to remove the nail polish, since it can affect the accuracy of the measurement. The sensor cannot be used inside or near the MRI. Avoid exposing the sensor to the intensive light source. In case of the long use, the place of application should be changed once every 4 hours to ensure the integrity of the patient's skin. In case of hyperactive and overly mobile patients, adhesive sensors should be used or applied with the use of adhesive tape.

## 4. Application

1. Connect the SpO<sub>2</sub> sensor to the specific pulse oximeter, turn the device on and check if it operates correctly, according to the manual.
2. Choose the place of application on the patient's body. Preferred application spots:
  - adults and children: index finger, alternatively thumb, middle finger, ring finger or a big toe
  - infants and neonates: big toe, alternatively palm, wrist, ankle or foot
3. Place the sensor on a patient's body and position infrared connection. Fastening should be lasting but should not cause pressure.
4. Regularly monitor the skin condition in the place of the SpO<sub>2</sub> sensor's application to ensure correct blood circulation and integrity of patient's skin.

After the SpO<sub>2</sub> sensor has been placed you have to ensure the cable is aligned in a way that prevents the patient from entangling or breaking the product. The SpO<sub>2</sub> adapters that simplify the application of the sensor and cable's attachments such as bands, handles or clips are available in manufacturer offer.

The safe and proper use of the SpO<sub>2</sub> sensor and its connection cable requires systematic checkups.

The checkups must be conducted at least couple of times a month depending on the frequency of use and cable's disinfection.

Visually examine the front of the isolation, connection joints and pins etc. Examine the mechanical integrity of the joint. The sensor should not be used if there seems to be any mechanical or electrical fault.

## 5. Suggestions for cleaning, disinfection, sterilization and instructions for use and storage

Cleanse the sensor and the connection cable with soap and warm water. Gently cleanse the sensor and the cable with the use of soft cloth moistened with water with soap or 70% isopropyl alcohol. Be careful not to damage the surface emitting the red and infrared radiation. The sensor and the cable should dry off before use. Do not use corrosive substances or chemicals other than 70% isopropyl alcohol. Do not soak or submerge the sensor. Excessive strength exercised during cleaning the sensor might cause the internal wires to tear. Those types of damage are not covered by the guarantee. Avoid any contact with strong dissolvents: aromatic, chlorinated, ketone, ether or ester. Do not allow the cleaners to have contact with the sensor's plug.

### Disinfection

The sensor should be disinfected with 70% isopropyl alcohol. Do not allow the alcohol to have contact with the plug from the joint side.

### Sterilization

The preferred method of sterilization is with ethylene monoxide (procedure according to gas sterilizer's manufacturers). After sterilization and before use the sensors must be precisely and safely aired. The other methods of sterilization might permanently damage the sensor.

## 6. Packaging and storage

The SpO2 sensors are packed individually. The sensors should be stored in the original pack according to the storage conditions so that its operating time could be as long as possible.

Storage conditions:

- temperature: from -25°C to +55°C
- relative humidity: from 10% to 85%, do not allow the steam to condense

The sensors and their cables, when not in use, should be loosely rolled up, without strong bends.

## 7. The conditions of the surroundings during operating

- temperature from +5°C to +45°C
- relative humidity from 30% to 75%
- atmospheric pressure 700 ÷ 1060 hPa

## 8. Warnings <sup>Δ</sup>

The SpO2 sensors are designed to be used with specific types of pulse oximeters. Before the application, the user is responsible for checking the compatibility of the sensor and the pulse oximeter. The incompatibility of those elements might reduce the accuracy and efficacy of the measurement and in an extreme case damage the product. The accuracy of the sensor should not be checked with functional tester or oximeter's simulator. The manual of a given device and its equipment should be checked to ensure their compatibility.

The use of the intravascular dye may distort the accuracy of the measurement. The incorrect application of the SpO2 sensor may distort the accuracy of the measurement as well. The nail polish and/ or tips should be removed before the application of the sensor since they may lead to inaccurate measurement. The SpO2 sensors should not be applied onto energized patient since their motion negatively affects measurement.

Do not use the SpO2 sensor during MRI or CT tests. The conducted electricity might cause burns, and MRI and CT devices might cause imprecise measurements. Also the sensor might affect the work of MRI and CT.

Strong electromagnetic field emitted by the electro-surgical equipment might affect the proper functioning of the SpO2 sensors.

Strong sources of visible and infrared light such as direct sunlight, light from the devices measuring the level of bilirubin or infrared heat radiation might affect the functioning of the sensor and as a result distort the readings. In case either of the situations occur, cover the spot of the sensor's application with opaque, non-transparent material.

Do not use other SpO2 bands or adapters than indicated by the manufacturer when placing the sensor or its cable on patient's body.

The additional pressure on the limb may result in alteration of venous pressure that may distort the findings.

The allergic reactions have not been observed when placing the sensor on the human skin, however such possibility cannot be entirely excluded, especially when it comes to allergically susceptible people.

In some cases there may be a need to regularly change the place of monitoring area, because the people skin's tolerance on the sensor might differ greatly.

Mobile radio-communication devices might influence the readings of the pulse oximeter. When the monitor does not display the reliable pulse measurement it may be caused by improper place of the sensor's application, or the monitored area is too thick, too thin or with deep pigmentation, which disables the normal light transmission. Than the place of the sensor's application must be changed.

The sensor that appears to have a mechanical or electrical defect, or its package is damaged must not be used.

Inundated sensors must not be used as well.

Do not alter or modify the sensor in any way. Such actions may influence the preceisness of readings.

The used sensor should be disposed of in accordance with local requirements.

## **8. Guarantee**

The guarantee lasts for one year, beginning at the day of the purchase.

Guarantee does not cover the improper use and modification of the sensors.

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The SpO<sub>2</sub> sensor generates, uses and may emit electromagnetic waves of radio frequency (RF). If the SpO<sub>2</sub> sensor is not used according to the directions specified in this instruction manual, it may cause electromagnetic interference. The user who uses the SpO<sub>2</sub> sensor should ensure that the sensor is used in an appropriate environment.

The SpO<sub>2</sub> sensors have been tested and recognized as compliant with restrictions specified in the EN 60601-1-2 standard for medical equipment. These restrictions can ensure suitable protection against electromagnetic interference (EMC), if the device is used in an environment which is suitable for the device and described in this instruction manual.

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The SpO2 sensor is intended for use in an electromagnetic environment which is described below. The buyer or user of the SpO2 sensor should ensure using the device in such environment.

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Note 1: In case of frequencies of 80 MHz and 800 MHz, higher frequency range should be used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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a. It is not possible to theoretically calculate precisely the electromagnetic field intensities generated by stationary transmitters, such as base stations for radio telephony (mobile/wireless), by land-based movable transmitters, amateur transmitters, FM and AM radio broadcasting stations and TV broadcasting stations. To define a state of an electromagnetic environment resulting from the operation of the stationary radio frequency transmitters, it is necessary to make measurements on site. If measured electromagnetic field intensity at the location where the SpO2 sensor is used exceeds the admissible level, you should monitor whether the SpO2 sensor is operating properly. If improper operation of the SpO2 sensor is observed, further actions may be necessary, such as turning or moving the device.

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The SpO2 sensor is intended for using in electromagnetic environments characterized by controlled interference resulting from radio frequency emission. The buyer or user of the SpO2 sensor may prevent electromagnetic interference by maintaining the minimum distance between portable and movable communication devices using radio frequency (transmitters), and the SpO2 sensor, according to the following recommendations, depending on the maximum output power of the communication equipment.

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Rated output power of the transmitter (W) Maximale

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For transmitters characterized by maximum rated output power which has not been mentioned above, a recommended separating distance in meters (m) can be determined based on the equation suitable for the transmitter frequency, where P is a maximum rated output power of the transmitter in watts (W), specified by the manufacturer of the transmitter.

Note 1: For frequencies of 80 MHz and 800 MHz, the separating distance for higher frequency range is used.

Note 2: These recommendations may not be applicable in all situations. Propagation of electromagnetic waves depends on absorption and reflections from buildings, objects and people.

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