SORIMEX

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Limb clamp electrodes

Extremitäten-Klammerelektroden

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EN

Instructions for use and storage of limb clamp electrodes

These instructions are for reusable EKK type ECG electrodes.

1. General characteristics

EKK limb electrodes are intended to receive bioelectrical signals from the human body, in order to transfer them to the ECG device.

Limb electrodes are reusable products.

2. Electrode characteristics

A limb electrode consists of:

- a clasp, intended to be put on the patient's limb (an arm or leg),
- an electroconductive part, in contact with the patient's skin,
- a fitted connection to be connected with ECG equipment.

The limb electrode is connected to the ECG equipment using an ECG cable fitted with a Ø 4mm diameter banana plug, type H latch, or type M clamp.

3. Guidelines of use

When using ECG limb electrodes, the following guidelines principle should be followed:

- the skin at the point of electrode placement should not be covered with an ointment, lotion or other substances,
- in places of electrode's contact with patient's skin, the skin surface should be cleaned with a swab moistened with 70% ethyl alcohol,
- electrodes are made in the following colours: black, red, yellow, green. Placement of the electrode of a specific colour on the appropriate limb should be done according to generally accepted principles or guidelines included in the ECG equipment instructions of use.
- do not open the electrode beyond its structural limit. If the patient's limb is of greater size than maximum opening of electrode jaws, use electrodes of a different structure or size. Forced placement of the electrode may lead to the cracking of the elastic element. Damage of this type is not covered by the guarantee.
- the connection cable should be arranged so as to prevent accidental movement of the electrode from the point of placement,
- in order to connect the electrode to ECG equipment, use cables structurally compatible with a given type of ECG equipment,
- wires connected to the electrodes should be protected against accidental contact with other conductive parts, as well as earthed elements.

4. Application

- a) Choose the point of application of the limb electrodes on the patient's body. Remove lotion, gel, ointments, dirt and other non-conducting layers from the chosen points of application. If necessary, degrease the patient's skin using water and soap or other appropriate washing/disinfecting products, and then carefully dry the skin before placement of the electrode. Do not apply electrodes on moist, lubricated or oily skin, for this would make it impossible to establish a proper connection and signal transmission.
- b) Apply the gel intended for ECG examinations on the patient's skin at the point of contact with the electroconductive part of the ECG electrode, using enough to ensure contact with the whole metal surface of the electrode.
- c) Connect the ECG cable connectors to the clasp electrodes.
- d) Place the electrodes on the patient's prepared skin.

5. Recommended cleaning, disinfection, and sterilization

These steps should be taken every time before using the electrodes for the examination of the next

patient. Before starting the cleaning or disinfection, the clasp electrodes should be disconnected from the ECG cable which they were connected to.

Cleaning

Wash the electrode with warm soapy water, delicately clean using a soft cloth. If necessary, use mild washing agents which will not react with the silver-covered electrode layer. The electrode should be dry before use. For washing, do not use agents containing abrasive materials, for example those intended to clean metals etc. When cleaning, do not soak or immerse the electrodes and do not wash under running water, do not use excessive strength or tools which may damage the electrode.

Disinfection

Clasp electrodes should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry.

Warning: do not use ultrasonic washers.

Sterilization

The product is not intended to be sterilized by any available methods.

6. Packaging and storing

Electrodes should be stored in their original packaging, in rooms at a temperature from 0°C to 50°C and air with relative humidity from 20% to 80%. Protect against sunlight, moisture and dirt.

Electrodes in use should be stored in conditions that have no negative effect on silver (in particular, in conditions that do not contain sulfur compounds).

7. Warnings

Electrodes are delivered to the client in a non-sterile state.

Before initiation of the ECG examination, it should be checked that the ECG device is functioning correctly. Use only together with compatible patient cables and ECG adapters.

Wires connected to the electrodes should be protected against accidental contact with other conductive parts, as well as earthed elements.

Electrodes should be connected to the ECG apparatus only by qualified medical personnel.

An ECG electrode is a product which has direct contact with the patient's skin; therefore allergic reactions may occur in patients prone to allergies. Occurrence of allergic reactions, or other skin-related reactions, requires immediate doctor's advice.

During the application of ECG electrodes, remember to apply sanitary procedures. Recommended methods of patient skin preparation before the application of the electrode: wipe the skin with several gentle movements with a gauze or swab moistened with a washing or disinfecting agent without injuring the epidermis, then dry the skin. In the case of damage to the epidermis during the procedure of preparation of the skin for electrode application, another point of application should be chosen.

Do not apply electroconductive gel and electrodes on injured or irritated skin.

Immediately after positioning and/or moving the patient, check the connection of electrodes with the ECG cables/wires and the ECG device is showing the appropriate reading.

Do not use the electrode if it is wet, found in damaged packaging, or for which the expiry date has passed. Excessive patient movement negatively influences the readings. Therefore, in the case of excitable patients, disruptions in ECG signals reception are possible.

Do not use ECG measurement electrodes immediately after defibrillation, for the time of signal recovery may be longer than 10 seconds.

Do not use ECG electrodes during Magnetic Resonance (MRI) examination or Computed Tomography (CT). Conducted electric current may inflict burns, and MRI and CT devices may cause inaccurate readings. An ECG electrode may influence MRI and CT images as well.

Under no circumstances should there be any contact with external electric current!

Disposal of used electrode should be done in accordance with legal requirements.

8. Shelf life

Prolonged and correct functioning of the electrode mainly depends on proper handling of the Ag/AgCl layer located on the electrode part adhering to the patient's body. The layer is applied in a galvanic process and requires careful handling. Even improper cleaning, e.g. with agents containing abrasive

substances, may damage this layer.

An electrode with damaged Ag/AgCl layer may receive electrical signals improperly and be a cause of improper ECG readings.

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