PL Elektrody EKG

EN ECG electrodes

DE EKG-Elektroden

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SK EKG elektródy

CS Elektrody EKG

RO Electrozi ECG

ET EKG elektroodid

LV EKG elektrodi

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EN

Instruction for disposable ECG electrodes, EK-S type.

1. Intended use and operation of the device

ECG electrodes are intended for diagnostic tests and monitoring of the patient's heart functions.

Electrodes are placed on the skin to transmit an electrical signal from the surface of the body via ECG cables and leadwires to the ECG apparatus or cardiac monitor, which convert the signal into electrocardiogram or vector cardiogram.

2. Patient target groups

Electrodes are dedicated to use by qualified medical personnel in infants, children and adults, as suggested by the manufacturer in point 12.

3. Contraindications

There are no known contraindications for use of the device.

4. Possible side effects

ECG electrode is a product that comes into direct contact with the patient's skin. There is a risk of irritation and allergic reactions in people with extremely sensitive skin.

The occurrence of skin reactions requires medical advice.

5. Warnings

- 5.1. Electrodes are disposable medical devices. Do not use again. Re-use may cause patients infection.
- 5.2. Use electrodes only on intact skin. Do not stick the electrode on wounded, infected or lesioned skin.
- 5.3. Do not use ECG electrodes during magnetic resonance imaging (MRI). Leaving electrodes on the skin may cause burns.
- 5.4. When performing electrosurgical procedures, place the electrodes as far away from the treatment site as possible to minimize the risk of burns to the patient's skin.

6. Cautions

- 6.1. Do not stick the electrodes on moist, lubricated skin due to the risk of the electrode coming loose.
- 6.2. Rough removal of the electrode may cause skin irritation.
- 6.3. When replacing an electrode with a new one, change the place of application of the next electrode to avoid skin irritation.
- 6.4. For disposable ECG electrodes, clamps and snaps are appropriate types of ECG cables/leadwires terminations. In case of cables/leadwires terminated with a banana plug, an adequate ECG adapter must be used.
- 6.5. Do not use electrodes that are wet or in damaged packaging.

- 6.6. Do not use electrodes which are expired.
- 6.7. Do not open the packaging when there is no necessity of using electrodes for the tests because of the risk of the gel drying out.
- 6.8. The use of alcohol to prepare the skin for testing is not recommended due to its skin impedance increasing effect. However, if alcohol is used, make sure it has completely evaporated from the skin before applying the electrode.
- 6.9. Excessive movement of the patient adversely affects the quality of the ECG signal.
- 6.10. In case the electrode comes off of the skin, replace it with the new one.
- 6.11. Applying electrodes on excessively hairy skin may be unable for proper ECG testing.
- 6.12. Do not use abrasive agents on infants, children and people with particularly sensitive skin.

7. Directions for use

- 7.1. Electrodes marked with „RTG” are intended for use in X-ray examination.
- 7.2. To connect the electrode with ECG apparatus, use ECG cables and leadwires compatible with equipment used.
- 7.3. Skin at the electrode site should be clean, dry and devoid from excessive hair.
- 7.4. If necessary, the skin should be shaved, washed with soap and water and dried thoroughly.
- 7.5. The oblong hole or cut in some types of electrodes is intended for attaching the leadwire.
- 7.6. Do not change the position of the electrodes once attached.
- 7.7. After the examination has ended, electrodes should be disposed according to the local regulations.
- 7.8. If a small amount of gel is left on the skin after removing the electrode, remove it with paper towel or wash it off with water.

8. Application

- 8.1. Open the pouch along the dotted line indicated on the packaging.
- 8.2. Determine the application sites on the patient's body. If necessary, prepare the patient's skin in accordance with point 7.3. Before sticking the electrode, always make sure that the application area is dry.
- 8.3. Connect the leadwire with electrode placed on the transparent foil.
- 8.4. After removing the electrodes from a transparent foil, apply to previously prepared skin, taking into account the information contained in points 5, 6 and 7. It is recommended to remove from the sticky side of the electrode the white paper element ("fingerlift"), which enables removing electrode from the foil. Apply the electrode to the skin and gently press its edges.

Fig. 1. Way of leadwire and electrode connection

8.5. For Holter electrodes:

8.5.1. Without removing the electrode from the transparent foil, pull the leadwire through the hole as shown on fig. 2, starting from the sticky side towards the outer side of the electrode.

8.5.2. Connect the leadwire located in the hole to the electrode. Orient the snap/clamp with the leadwire upwards as shown in fig. 2. Leave approx. 3-4 centimeters of loose wire (above the hole) to prevent it from disconnecting from the electrode.

8.5.3. Make sure to stick the electrodes with the hole for attaching the leadwires facing downwards (fig. 2).


Fig. 2. Way of leadwire an Holter electrode connection

9. Electrode removal

After the examination, lift the electrode and gently, slowly detach it while holding the patient's skin with your free hand. Avoid peeling off at right angles to minimize the risk of skin irritation.


10. Storage conditions

10.1. Electrodes should be stored in their original packaging, at a temperature between 5°C and 30°C. Protect from direct sunlight, moisture and dirt.

10.2. The expiry date of the product is given on the packaging and applies to originally sealed products, stored in accordance with the manufacturer's recommendations.

10.3. Leave unused electrodes in their original packaging and use within the next 7 days. Secure the packaging by folding its edges.

11. Explanation of symbols

| | | | | | |
|--|---|---|---|--|---|
|  | Suggested group of patients: adults |  | Indicates the date |  | Protect against sunlight |
|  | Suggested group of patients: pediatric. Applies to children from infancy (2 years) to 18 years of age, depending on the patient's weight. |  | Production batch number |  | Keep dry |
|  | |  | Product UDI code |  | Temperature limit |
|  | Liquid gel |  | Manufacturer |  | Do not re-use |
|  | Solid gel |  | CE mark (conformance to the Regulation (EU) 2017/745) |  | Caution |
|  | Catalogue number |  | Medical device |  | Product has an Ag/AgCl sensor |
|  | Product index |  | Read the instruction for use |  | Product does not contain natural latex |
|  | Pieces per packaging |  | Do not use if the package is damaged |  | Product does not contain polyvinyl chloride |
|  | Importer |  | Authorized representative in the Switzerland | | |

12. Suggested patient groups and intended use of ECG electrodes type EK-S

| REF | Patient group | Intended use |
|-----------------|---------------|--------------------------------------|
| EK-S 22 PSG | pediatric | resting examination |
| EK-S 25 PSG | pediatric | short-term monitoring |
| EK-S 30 PSG RTG | pediatric | short-term monitoring |
| EK-S 30 PSG | pediatric | resting examination |
| EK-S 35 PSG ST | pediatric | stress tests |
| EK-S 36 P | pediatric | long-term monitoring |
| EK-S 39 PSG | pediatric | Holter |
| EK-S 40 PSG | pediatric | long-term monitoring |
| EK-S 42 PSG | pediatric | long-term monitoring |
| EK-S 44 PSG | for adults | resting examination |
| EK-S 45 PSG RTG | for adults | short-term monitoring |
| EK-S 45 PSG | for adults | short-term monitoring |
| EK-S 48 PSG | for adults | stress tests |
| EK-S 50 NWSG | for adults | Holter |
| EK-S 50 P | for adults | long-term monitoring |
| EK-S 50 PSG | for adults | long-term monitoring |
| EK-S 55 P | for adults | short-term monitoring / stress tests |
| EK-S 55 PSG | for adults | long-term monitoring |
| EK-S 56 P | for adults | resting examination |
| EK-S 60 P T | for adults | Holter |
| EK-S 60 P | for adults | Holter |
| EK-S 60 PSG XC | for adults | stress tests |
| EK-S 60 PSG | for adults | Holter |
| EK-S 61 P | for adults | long-term monitoring |
| EK-S 61 PSG | for adults | long-term monitoring |
| EK-S 62 PSG XC | for adults | stress tests |

13. Product incident reporting

Any serious incident related to the device should be reported to the manufacturer and to the authority competent for the place of residence of the user or patient.

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