# SORIMEX

## Kable SpO2 SpO2 cables

SORIMEX sp. z o.o. sp. k.  
ul. Równinna 25, 87-100 Toruń, Poland  
tel.: +48 56 657 77 20, fax: +48 56 657 77 21  
e-mail: [biuro@sorimex.pl](mailto:biuro@sorimex.pl) ♦ [www.sorimex.eu](http://www.sorimex.eu)

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EN

## Instruction of use and storage SpO2 cables

### 1. General description

The SpO2 cable is intended for the electrical connection of a saturation sensor with a device for measuring oxygen saturation (pulse oximetry). SpO2 cables are classified as Class I and delivered as non-sterile.

### 2. Characteristics

The SpO2 cable consists of a suitable multi-core cable and its connector is adapted to the saturation sensor and the device to which it is to cooperate. The user should make sure that the cable is suitable for the specific type of saturation sensor and device .

Cables are intended for use with undamaged patient's skin.

### 3. Principles of use

When using SpO2 cables, the following principles should be observed:

- make sure that SpO2 cable is suitable for the device to work with,
- cables should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the cable,
- place the cable so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin) when connecting and disconnecting the cable,
- protect the cable connectors against suffusion.

### 4. Installation

Connect the SpO2 cable from one side to the SpO2 sensor and the other side to the device. Turn on the unit and check that the monitor operates properly according to the monitor's instruction manual. Pay attention to the proper selection of the SpO2 cable for the appropriate device and SpO2 sensor.

## **5. Recommended cleaning, disinfection and sterilisation**

These activities should be realized every time before use of the device and cable before the examination of the next patient. Before starting the cleaning, disinfection, or sterilisation, the cable should be disconnected from the device with which it cooperates.

### Cleaning

Wash the cable with warm soap water, delicately clean using soft cloth wetted with soap water (mild detergent) or 70% isopropyl alcohol. The cable should dry before use.

Do not use corrosive substances or other chemical preparations.

When cleaning, do not soak or immerse the cables, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to cable damage. Continuous bending of the conductors during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.

Avoid contact with strong solvents: aromatic, chloric, ketone, etheral, and ester. Do not allow for contact of cleaners with tplugins and pins on either side of the cable.

### Disinfection

Cables should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry. Do not allow for contact of cleaners with tplugins and pins on either side of the cable.

### Sterilization

The recommended method consists of sterilization using ethylene oxide. After sterilization and before use, cables should be accurately degassed in device intended for that purpose.

Sterilization with other methods is not recommended for it may permanently damage the cable.

## **6. Packaging and storing**

SpO2 cables are packed individually. In order to maximize the product's working life, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:

- temperature: -25°C to + 55°C,
- protect against sunlight, moisture and dirt.

Cables in use should be stored in a dry place.

When not used, SpO2 cables should be loosely coiled, without sharp bends.

## **7. Warnings**

Failure to observe the above rules may result in a break in the electrical circuit in the transmission cable or cause a short circuit. These symptoms are easy to identify by operating the camera. The damaged cable must be replaced and replaced with a working one.

Utilisation of used cable should be done in accordance with binding provisions.

## **8. Guarantee**

The guarantee period covers 1 year from the date of sale.

The guarantee does not include defects resulting from improper use and storage.

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